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EN - INDEX

	a The Loss development	
7.	ALERTS 1	0
6.	FUNCTIONS AND NAVIGATION	8
c l		~
5.	GENERAL INFO	7
1.	INTRODUCTION	7

1. Introduction

Thank you for choosing our equipment. This manual describes the controller main functions to provide the best cooling conditions, in compliance with the applicable European directives and standards, with the aim of providing a safe and practical use that will meet your needs and expectations.

5. General info

On / Off

With the equipment connected to the power mains and in stad by mode, press the On / Off button to turn it on. To log in with your credentials, press the marked area and enter your password.

♪

To turn off the equipment, press the On / Off button in the main menu. Do not disconnect the equipment from the power supply as this will be interpreted as a power failure, generating

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6. Functions and navigation

Operating temperature

Users allowed to change the operating temperature can do so from this menu. To change press the keys and until you reach the desired temperature. Press store and return to main menu.

Manual defrost

Active a defrost cycle from this menu is the device conditions allows it.

+ -

Battery status

Check battery status from this menu. If there is a problem with the battery performance it will be displayed "Battery problems". If there is a communication problem it will be displayed "information unavailable". In the main menu consult the battery charge status.

Door lock

Only available for devices supplied with electronic locks. This buttons allows to lock the door.

Γġ

Glass door frame heater

The glass door is equipped with a frame heater for anti-fogging. To activate this function proceed as follows according to your model.

Acura models:

- 1. Access the Settings menu;
- Press Set Parameters and entre password 743;
- Scroll down and select parameter u06c;
- 4. Set **u06c = 6**;

5. To save, press the set parameter again; The equipment will adapt its operation and whenever necessary activate the resistance of the Exige models:

- 1. Open the frontal cover to access the equipment cooling unit;
- 2. Turn On the door frame heater button.

The equipment will adapt its operation and whenever necessary activate the resistance of the

glass door.

- History
- This menu displays the equipment historical operation data of:
 - ✓ Alerts;
 - ✓ Defrost;
 - ✓ Door openings.

To restart the count press the button .

To delete the data press the button and insert the password **4-5-9-8**.

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6.1.1. Alerts history

Consult the historical data from the past 30 alerts, as follows:

- Total number;
- Type;
- ✓ Start and finish date / time;
- ✓ Max. Temperature during alert.

. ۲

6.1.2. Defrost history

Consult the historical data from the past 30 defrosts, as follows:

- ✓ Total number;
- Type;
- \checkmark Start and finish date / time.

6.1.3. Door open history

Consult the historical data from the past 30 door openings, as follows:

- ✓ Total number;
- ✓ Critical number;
- ✓ Total time.

6.1.4. Equipment history

Consult a report of the equipment performance, as follows:

- ✓ Total hours of equipment On (Hr);
- ✓ Total hours of compressor On (Hr);
- ✓ Average time in minutes of compressor On (min.);
- ✓ Average time in minutes of compressor Off (min.);
- ✓ Average time percentage of equipment On (%);
- ✓ Average time percentage of compressor On in the last hour (%).

6.1.5. History settings

Manage and set the historical data you would like to consult, from the following:

- ✓ Temperature probe;
- ✓ Temperature product probe;
- ✓ Temperature evaporator probe;
- ✓ Temperature condenser probe;
- ✓ Temperature reference probe;
- ✓ Temperature regulation probe;
- ✓ Operating temperature;
- ✓ Equipment status;
- ✓ User log in;
- ✓ User actions;
- ✓ Defrost;
- ✓ Door open;
- ✓ Alerts.

Temperature alerts

Manage pre-set alerts. Browse with the \square and \square buttons to check and set minimum and maximum threshold for the following alerts:

- ✓ Operating temperature;
- ✓ Product probe temperature;
- ✓ Pre-alert operating temperature;
- ✓ Pre-alert probe temperature.

User management

Manage and set user access levels. Set up to 9 user access permissions (1 Admin + 8 users). To manage users you must be log in as a user Admin, by default the password to log in as a user Admin is **0-0-0-0**. The user Admin management permissions are:

- ✓ Activate / deactivate users;
- ✓ Manage user access permissions;
- ✓ Manage and set log in passwords per user.

To change the user Admin default password, log in and access the user management menu and press the code area to edit.

ENT

Language

Language settings menu, available languages are: English, Italian, Spanish, German, French, and Portuguese.

Ŋ

Service

Consult and manage the following:

- ✓ Set date and time;
- Consult inputs and outputs;
- ✓ Consult alerts list;
- ✓ Set equipment name;
- ✓ Manage welcome page;
- ✓ Manage advanced functions.

21

6.1.6. Advanced functions

Perform probes calibration and do performance tests.

Settings

Set and manage equipment parameters and restore default parameters.

7. Alerts

When an alert is active the buzzer will sound and this icon will be displayed on the screen. There are alerts that can be automatically silenced, if the problem is corrected within an accepted timeframe. Other alerts need to be manually silenced.

Alert ID	Description
CABINET PROBE AL.	Control probe alert. Motives: probe disconnected and/or faulty. Actions: replace probe; contact service maintenance.
EVAPORATOR PROBE AL.	Defrost probe alert. Motives: probe disconnected and/or faulty. Actions: replace probe; contact service maintenance.
CONDENSER PROBE AL.	Condenser probe alert. Motives: probe disconnected and/or faulty. Actions: replace probe; contact service maintenance.
PRODUCT PROBE AL.	Product probe alert. Motives: probe disconnected and/or faulty. Actions: replace probe; contact service maintenance.
REFERENCE PROBE AL.	Reference probe alert. Motives: probe disconnected and/or faulty. Actions: replace probe; contact service maintenance.
DOOR OPEN	Open door alert. Motives: door opened for a long period of time; door is not closing properly; reed sensor faulty. Actions: inspect the problem; replace reed sensor; contact service maintenance.
TIMEOUT DEFROST	Alert of defrost time exceeded. Actions: perform manual defrost; turn off the equipment and clean the excess of ice; contact service maintenance.
HIGH TEMP. CABINET	High internal temperature alert. Actions: inspect door for damage; avoid having the door open for long periods of time; contact service maintenance.
PRE-ALARM HIGH T. CABINET	High internal temperature pre-alert. Actions: inspect door for damage; avoid having the door open for long periods of time; contact service maintenance.
LOW TEMP. CABINET	Low internal temperature alert. Motives: unconfigured parameters; faulty microprocessor. Actions: configure parameters; replace the microprocessor; contact service maintenance.
HIGH TEMP. PRODUCT	High product probe temperature alert. Motives: probe disconnected and/or faulty; unconfigured parameters; communication with <i>Backup</i> module faulty; contact service maintenance.
PRE-ALARM HIGH T. PRODUCT	High product probe temperature pre-alert.
PRE-ALARM LOW T. PRODUCT	Low product probe temperature pre-alert.
Low Temp. Product	Low product probe temperature alert. Motives: unconfigured parameters; faulty microprocessor. Actions: configure parameters; replace the microprocessor; contact service maintenance.

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