

Eckert & Ziegler Co-57 FLOOD SOURCE APPLICATION CHART

Camera Mfr.	Camera Model	SPECT or PLANAR	Source Type	Head Size	Recommended Nominal Acvity	Perflexion™	Feather Lite™	NES	C-Thru
DIGIRAD	2020tc	SPECT	Square	8"x 8"	10mCi	PF09R-057-10M	MED3722	N/A	CTRFI0012
	Spect PAK	SPECT	Square	9"x9"	10mCi	PF09R-057-10M	MED3722	N/A	CTRFI0012
	Cardius Series	SPECT	Square	8" x 6"	10mCi	N/A	MED3722	N/A	CTRFI0012
Dilon (Preferred OEM SPEC)	6800	SPECT	Rectangular	8" x 6"	5mCi	PF16R-057-5M	MED3715	NES8505	CTRFI0012
ELSCINT	SP4	SPECT	Circular	15.75"	10mCi MAX	PF18C-057-10M	MED3706	NES392	CTRFI0050
	SP6	SPECT	Rectangular	21.25"x15.75"	10mCi MAX	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Helix	SPECT	Rectangular	21.25"x15.75"	10mCi MAX	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Cardial	SPECT	Rectangular	15.75" x 10"	10mCi MAX	PF16R-057-10M	MED3713	NES8510	CTRFI0017
	Varicam	SPECT	Rectangular	21.25"x15.75"	10mCi MAX	PF24R-057-10M	MED3709	NES8400	CTRFI0000
GENERAL ELECTRIC MEDICAL SYSTEMS	300 AM	PLANAR	Circular	10"	3mCi	N/A	MED3700	NES298	NA
	Maxi 2	PLANAR	Circular	15"	5mCi	PF18C-057-10M	MED3706	NES391	CTRFI0051
	Maxi 37	PLANAR	Circular	15"	5mCi	PF18C-057-10M	MED3706	NES391	CTRFI0050
	5000 SERIES	SPECT	Circular	24"	10mCi	N/A	MED3708	NES8012	NA
	Discovery VH	SPECT	Rectangular	20" x 14"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Discovery 530	SPECT	Rectangular	14.6" x 7.4"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	Hawkeye	SPECT	Rectangular	20" x 14"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	DSTI/Dsi	SPECT	Rectangular	12.99" x 14"	10mCi	N/A	MED3713	NES8480	NA
	Maxxus	SPECT	Rectangular	21" x 16"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Starcam XRT	SPECT	Rectangular	21" x 16"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Starcam XCT	SPECT	Circular	15.4"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
	Starcam ACT	SPECT	Circular	15.4"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
	Optima	SPECT	Rectangular	14" x 9"	7.5mCi	PF16R-057-7.5M	MED3727	NES8510	CTRFI0017
	Millenium MPS	SPECT	Square	14" x 14"	10mCi	N/A	MED3725	NES8480	NA
	Millenium MPR, VG	SPECT	Rectangular	15.75" x 21.75"	10 mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Millenium MG	SPECT	Rectangular	14" x 20"	10 mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Neurocam	SPECT	Rectangular	7.9" x 6.7"	10mCi	PF09R-057-10M	MED3722	N/A	CTRFI0012
	Myosite	SPECT	Rectangular	20" x 14"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Infinia	SPECT	Rectangular	21.25" X 15.75"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
Ventri	SPECT	Rectangular	14.6" x 7.4"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017	
HITACHI	1024C	SPECT	Circular	13.5"	5mCi	N/A	MED3700	NES298	N/A
	1024 RDT	SPECT	Rectangular	19.7" x 14.2"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	1024 R	SPECT	Rectangular	19.7" x 14.2"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Neuro Spect	SPECT	Rectangular	8.7" x 6.7"	10mCi	PF09R-057-10M	MED3722	N/A	CTRFI0012
	150 250DSP	SPECT	Rectangular	20" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	260 DSP	SPECT	Rectangular	22" x 16"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
IS2	SR (Single Head Rect.)	SPECT	Rectangular	20" x 14.8"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Pulse CDC	SPECT	Rectangular	15" x 10.2"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	DCC (Digital cardiac camera)	SPECT	Rectangular	15" x 10.2"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	SC (Single Head Circular)	SPECT	Circular	15"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
MEDISO MEDICAL (NC Systems)	CardioSpect SC	SPECT	Rectangular	20.8" x 15.3"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	CardioSpect SR	SPECT	Rectangular	20.8" x 15.3"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	CardioSpect D90	SPECT	Rectangular	14.5" x 9"	10mCi	N/A	N/A	N/A	CTRFI3770A
	CardioSpect VMAX	SPECT	Rectangular	20.8" x 15.3"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	NeuroSpect Quad	SPECT	Rectangular	9" x 8"	10mCi	PF09R-057-10M	MED3722	N/A	CTRFI3770A
NUCLEAR CHICAGO	Phogamma	PLANAR	Circular	12"	5mCi	N/A	MED3700	NES297	N/A
	LEM	PLANAR	Circular	9.75"	3mCi	N/A	MED3700	NES297	N/A
	LFOV	PLANAR	Circular	15"	5mCi	PF18C-057-5M	MED3706	NES391	CTRFI0051
PARK MEDICAL	Isocam I & II	SPECT	Rectangular	22.3" x 16.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000

Camera Mfr.	Camera Model	SPECT or PLANAR	Source Type	Head Size	Recommended Nominal Acvity	Perflexion™	Feather Lite™	NES	C-Thru
PHILIPS formerly ADAC	ARC 3000	PLANAR	Circular	15"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
	Genesys	SPECT	Rectangular	20"x15"	10mCi MAX	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Brightview	SPECT	Rectangular	20.25" x 16"	10mCi MAX	PF24RW-057-10M	N/A	NES8610	CTRFI0021
	Argus	SPECT	Rectangular	20"x15"	10mCi MAX	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Solus, Cardial, or Vertex	SPECT	Rectangular	20"x15"	10mCi MAX	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Forte	SPECT	Rectangular	20"x15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Cirrus	SPECT	Circular	15"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
	Cardio MD	SPECT	Rectangular	9.2" x 15.4"	10mCi MAX	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	Cardio 60	SPECT	Rectangular	20" x 15"	10mCi MAX	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Precedence	SPECT/CT	Rectangular	20x15	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
Skylight	SPECT	Rectangular	20" x 14"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000	
PHILIPS formerly Marconi/Picker	Dynamo	PLANAR	Circular	10"	3mCi	N/A	MED3700	NES297	N/A
	411	PLANAR	Circular	11"	5mCi	N/A	MED3700	NES297	N/A
	412	PLANAR	Circular	12"	5mCi	N/A	MED3700	NES297	N/A
	415	PLANAR	Circular	15"	5mCi	PF18C-057-10M	MED3706	NES391	CTRFI0051
	Prism 1000	SPECT	Rectangular	20" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Prism 2000	SPECT	Rectangular	20" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Prism 3000	SPECT	Rectangular	15.7" x 9.4"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	SX300	SPECT	Square	14" x 14"	10mCi	N/A	MED3725	NES8480	CTRFI0000
Axis / Irix	SPECT	Rectangular	21" x 15.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000	
RAYTHEON		PLANAR	Circular	15"	5mCi	PF18C-057-10M	MED3706	NES391	CTRFI0051
SIEMENS MEDICAL SYSTEMS	Multispect 2	SPECT	Rectangular	21.25" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Multispect 3	SPECT	Rectangular	16" x 12"	10mCi	N/A	MED3713	NES8480	N/A
	3700, 7500 Orbiter Series	SPECT	Circular	15.25"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
	Body Scan	SPECT	Rectangular	23.5" x 15.75"	10 mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Diacam	SPECT	Rectangular	21.25" X 15"	10 mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	3700, 7500 Orbiter Series	SPECT	Rectangular	21.25" X 15"	10 mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	c.cam	SPECT	Rectangular	14" x 8.4"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	Symbia*	SPECT	Rectangular	21.23" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Symbia Auto QC Source	SPECT	Line Source	21.23" x 15"	10mCi	NA	N/A	HEGL-0133/PHI-0124(point source)	N/A
	e.cam	SPECT	Rectangular	21.25" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
SMV (SOPHA MEDICAL)	DST	SPECT	Rectangular	11.7" x 15.8"	10mCi	N/A	MED3713	NES8480	N/A
	DSX	SPECT	Rectangular	21.3" x 15.75"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	DS7	SPECT	Circular	15.75"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
	Bodytrac	SPECT	Rectangular	21.3" x 15.75"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	DSTXL	SPECT	Rectangular	21.3" x 15.75"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Vision FX Series	SPECT	Rectangular	20" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
TECHNICARE (G.E.)	Mobile 420	PLANAR	Circular	10"	3mCi	N/A	MED3700	NES297	N/A
	Mobile 420	PLANAR	Circular	12"	5mCi	N/A	MED3700	NES297	N/A
	Mobile 420	PLANAR	Circular	15"	5mCi	PF18C-057-5M	MED3700	NES391	CTRFI0051
	Omega 500	SPECT	Rectangular	20" x 14.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Gemini 700	SPECT	Rectangular	20" x 14.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Gemini 600	SPECT	Circular	15"	10mCi	PF18C-057-10M	MED3706	NES392	CTRFI0050
TOSHIBA AMERICA	GCA7100A & 7200	SPECT	Rectangular	21.5" X 19"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	GCA 9300	SPECT	Rectangular	16.1" x 8.3"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	GCA901, 901A	SPECT	Rectangular	20.1" x 14.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	GCA901WB	SPECT	Rectangular	20.1" x 14.6"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	GCA602, 602A	SPECT	Circular	13.75"	5mCi	N/A	MED3700	NES297	N/A
	GCA601A	SPECT	Circular	13.75"	5mCi	N/A	MED3700	NES297	N/A
	t.cam	SPECT	Rectangular	21.25" x 15"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
TRIONIX RES. LAB.	Monad	SPECT	Rectangular	24" x 14.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Biad 24/XLT24	SPECT	Rectangular	24"x 14.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Biad 20/XLT20	SPECT	Rectangular	20" x 14.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
	Triad XLT	SPECT	Rectangular	15.75" x 8.7"	10mCi	PF16R-057-10M	MED3727	NES8510	CTRFI0017
	Triad XLT 20	SPECT	Rectangular	20" x 14.5"	10mCi	PF24R-057-10M	MED3709	NES8400	CTRFI0000
UNION CARBIDE	Union Carbide	PLANAR	Circular	16"	5mCi	PF18C-057-5M	MED3706	NES391	CTRFI0051