

Supplements for Dehydrated Culture Media

Item #	Supplement & Description	Package*
7980	Buffered Listeria Enrichment Supplement (FDA) , 5 mL Acriflavine, 4.5 mg Nalidixic acid, 18 mg Cycloheximide, 22.5 (Each vial supplements 450 mL of Buffered Listeria Enrichment Broth Base)	10 vials/pkg lyophilized
7981	Campylobacter Supplement (CFP) , 5 mL Cefoperazone, 16 mg (Each vial supplements 500 mL of Campy Blood Free Selective Medium)	10 vials/pkg lyophilized
7982	Egg Yolk Emulsion , 100 mL Chicken egg yolk and egg, 50% solution (Each bottle supplements 2 L of Baird Parker Agar)	1 each liquid
7983	Egg Yolk Tellurite , 100 mL Chicken egg yolk and egg, 100% solution Potassium Tellurite, 0.21 g (Each bottle supplements 2 L of Baird Parker Agar)	1 each liquid
7984	Fraser Broth Supplement , 5 mL Ferric ammonium citrate, 250 mg (Each vial supplements 500 mL of Fraser Broth Base)	10 vials/pkg liquid
7985	Novobiocin Supplement , 10 mL Novobiocin, 20 mg (Each vial supplements 1 L of EC Medium, Modified, 1 L of Rappaport-Vassiliadis (MSRV) Medium Semisolid, Modified, or 1 L of Tryptic Soy Broth, Modified)	10 vials/pkg lyophilized
7986	Oxford Listeria Supplement , 5 mL Cycloheximide, 200 mg Colistin, 10 mg Acriflavine, 2.5 mg Cefotetan, 1 mg Phosphomycin, 5 mg (Each vial supplements 500 mL of Oxford Listeria Agar)	10 vials/pkg lyophilized
7987	PALCAM Supplement , 5 mL Polymyxin B, 5 mg Ceftazidime, 10 mg Acriflavine, 2.5 mg (Each vial supplements 500 mL of PALCAM Media)	10 vials/pkg lyophilized
7988	Rosolic Acid 1% , 10 mL (Each vial supplements 1 L of m-FC Agar and m-FC Broth)	10 vials/pkg lyophilized
7989	Tellurite Supplement (1%) Chapman , 10 mL Potassium Tellurite, 100 mg (Each vial supplements 1 L of Baird Parker Agar)	10 vials/pkg lyophilized
7990	XLT4 Supplement , 100 mL Tergitol 4 (27%) (4.6 mL will supplement 1 L of XLT4 Agar)	1 each liquid

*All supplements stored at 4–8° C.