


Acumedia®

Microbiological Matrices

E. coli

	<i>Nonselective Media</i>	<i>Purity Plate</i>	<i>Selective-Differential</i>	<i>E. coli O157:H7</i>
Brain Heart Infusion Agar (7115)	X			
Brilliant Green Bile Agar (7449)			X	
EC Medium with MUG (7361)			X	
EC Medium, Modified (7506)				X
Eosin Methylene Blue Agar (7134)			X	
Lauryl Sulfate Broth with MUG (7300)			X	
MacConkey Agar (7102)			X	
MacConkey Agar with Sorbitol (7320)				X
Nutrient Agar (7145)		X		
Tryptic Soy Agar (7100)	X	X		
Violet Red Bile Agar (7165)			X	
Violet Red Bile Agar with MUG (7359)			X	



MacConkey Agar with Sorbitol



Violet Red Bile Agar

Coliforms

	<i>Food Analysis</i>	<i>Environmental</i>
A-1 Medium (7601)		X
Brilliant Green Bile Agar (7449)	X	X
Brilliant Green Bile Broth 2% (7119)	X	X
Brilliant Green Bile Broth 2% with MUG (7344)	X	X
EC Medium (7206)	X	X
EC Medium, Modified (7506)	X	
EC Medium with MUG (7361)	X	X
Eosin Methylene Blue Agar, Levine (7103)	X	X
Lactose Broth (7141)	X	X
Lauryl Sulfate Broth (7142)		X
Lauryl Sulfate Broth w/MUG (7300)		X
MacConkey Agar (7102)	X	X
MacConkey Agar with Sorbitol (7320)	X	X
MacConkey Agar, CS (7391)	X	X
MacConkey Agar, Modified (7440)	X	
MacConkey Broth (7185)	X	X
Tergitol 7 Agar (7187)	X	X
Violet Red Bile Agar (7165)	X	X
Violet Red Bile Agar with MUG (7359)	X	X
Violet Red Bile Glucose Agar (7425)	X	