

Product Information and Prices



Elmasonic Accessories - Practicals

Sieve cleaning module SRH 4 / 200

- made of stainless steel and with integrated rotating device
- up to 4 test sieves [Ø up to 200 / 203 mm (8")] can be cleaned simultaneously in the ultrasonic bath with the SRH 4 / 200. The inclined position of the sieves in the SRH 4 / 200 ensures a better cleaning and rinsing effect
- for gentle cleaning of the test sieves without changing the mesh size or fabric tension. This ensures that the dimensional accuracy and calibration of the sieves is maintained 100 %



	300	450	500	900
Order number 230V	1042706			
Order number 115V	1044904			
Dimensions W / D / H (mm)	675 / 295 / 330			
Weight (kg)	6,8			
Suitable for	300	450	500	900
• Elmasonic Easy	1 x			
• Elmasonic Med	1 x			
• Elmasonic Select	1 x	1 x		2 x
• Elmasonic S	1 x		1 x	2 x
• Elmasonic P	1 x			

Dimension definition

