

Brilliant Performance. Clear Results.

Seeing Beyond Limits



NEO
IRIS®



High throughput full automation

Maximum productivity

Performance

NEO Iris® automates Immucor's patented Capture® technology with proven clinical performance and a comprehensive process control strategy, providing trustworthy test results.

Productivity


NEO Iris delivers the highest type-and-screen throughput on the market. True continuous access allows operators to add samples or resources anytime.




- A** Windows-based user friendly touchscreen interface
- B** Linear rack sample and reagent bay; load up to 224 samples at once
- C** Dual pipettors for processing multiple samples and plates simultaneously
- D** Microplate loading tower, washer module, incubator module, image analysis reader
- E** Centrifuge module
- F** Fluidics module is self-contained

Clear results Remarkable flexibility


Dynamic scheduler view
with real time to results



Flexibility
NEO Iris' dynamic scheduler,
with STAT priority handling offers
superior workflow management.



Configurable Reflex Testing for
optimized lab workflow



Results screen with improved
image resolution


Broad Test Menu

- ABO/Rh D Typing
- Weak D
- Donor Confirmation
- CMV
- IgG DAT
- IgG Crossmatch
- Antigen Screening
- Phenotyping (Rh, Kell, Jka, Jkb, Fya, Fyb, S, s, k)
- Antibody Identification (Primary Panel, D Pos Panel & D Neg Panel)
- Antibody Screening (2-cell, 3-cell, & Pooled)

Brilliant performance Exceptional support

Integration made easy. NEO Iris is already integrated into the existing Immucor complete solutions offering, combined with supporting software, ImmuLINK® & blud-directSM.

Data Management



Remote Diagnostics

blud_direct

Blood Bank Labs can't afford to be down. That's why there's blud_direct. An industry exclusive premium level of secure, technical support utilizing high-speed Internet connectivity enabling you to receive help from an Immucor Technical Service Representative directly on the instrument monitor. With an overwhelming success rate, downtime is minimized instantly.

IMMUCOR

For more information about the product and to know more about the regulatory status in your country, please contact your local Immucor Representative, imarketing@immucor.com or visit www.immucor.com